
memo

Release 0.1.0

Arnaud Gaudry

Feb 18, 2022

CONTENTS:

1	memo_ms package	1
1.1	Submodules	3
1.1.1	memo_ms.classes module	3
1.1.2	memo_ms.import_data module	4
1.1.3	memo_ms.visualization module	5
2	Description	7
3	To install it:	9
4	Examples	11
5	Indices and tables	13
	Python Module Index	15
	Index	17

MEMO_MS PACKAGE

```
class memo_ms.FeatureTable(path: str, software: str)
```

Bases: `object`

Create a FeatureTable dataclass object from a feature table

Args: path (str): Path to a feature table file (.csv) software (str): One of [mzmine, xcms, msdial, memo]: the software used for feature detection.

Returns: self.feature_table (DataFrame): A cleaned feature quantification table

```
__init__(path: str, software: str) → None
```

```
export_matrix(path, sep=',')
```

Export a given matrix

Args: path (str): path to export sep (str): separator

Returns: None

```
filter(samples_pattern, max_occurence=None)
```

Filter a feature table: remove samples matching samples_pattern AND remove features occurring in more than n = max_occurence samples matched by samples_pattern

Args: samples_pattern (str): the str pattern to match in samples to filter max_occurence (int): maximal number of occurrence allowed in matched samples before removing a feature/word

Returns: self.filtered_feature_table (DataFrame): A filtered feature table

```
class memo_ms.MemoMatrix
```

Bases: `object`

Create an empty MemoMatrix dataclass object

```
__init__() → None
```

```
export_matrix(path, sep=',')
```

Export a given matrix

Args: path (str): path to export sep (str): separator

Returns: None

```
filter(samples_pattern, max_occurence=None)
```

Filter a MEMO matrix: remove samples matching samples_pattern AND remove features occurring in more than n = max_occurence samples matched by samples_pattern

Args: samples_pattern (str): the str pattern to match in samples to filter max_occurence (int): maximal number of occurrence allowed in matched samples before removing a feature/word

Returns: self.memo_matrix (DataFrame): A filtered feature table matrix

memo_from_aligned_samples(featuretable, spectradocuments) → pandas.core.frame.DataFrame

Use a featuretable and a spectradocuments to generate a MEMO matrix. Returns a pd.DataFrame MEMO matrix.

Args: featuretable (FeatureTable): a FeatureTable dataclass object spectradocuments (SpectraDocuments): a SpectraDocuments dataclass object

Returns: self.memo_matrix (DataFrame): A MEMO matrix

memo_from_unaligned_samples(path_to_samples_dir, min_relative_intensity=0.01, max_relative_intensity=1.0, min_peaks_required=10, losses_from=10, losses_to=200, n_decimals=2)

Generate a Memo matrix from a list of individual .mgf files

Args: path_to_samples_dir (str): Path to the directory where individual .mgf files are gathered
min_relative_intensity (float): Minimal relative intensity to keep a peak
max_relative_intensity (float): Maximal relative intensity to keep a peak
min_peaks_required (int): Minimum number of peaks to keep a spectrum
losses_from (int): minimal m/z value for losses
losses_to (int): maximal m/z value for losses
n_decimals (int): number of decimal when translating peaks/losses into words

Returns: self.memo_matrix (DataFrame): A MEMO matrix

merge_memo(memomatrix_2, drop_not_in_common=False)

Merge 2 MEMO matrix

Args: memomatrix_2 (MemoContainer): MemoMatrix dataclass object containing the 2nd MEMO matrix to merge
drop_not_in_common (bool): Drop peaks/losses not in common

Returns: MemoContainer (MemoContainer): A MemoMatrix dataclass object containing the merged MEMO matrix

class memo_ms.SpectraDocuments(path: str, min_relative_intensity: float = 0.01, max_relative_intensity: float = 1.0, min_peaks_required: int = 10, losses_from: int = 10, losses_to: int = 200, n_decimals: int = 2)

Bases: `object`

Create a SpectraDocuments dataclass object containing spectra documents and metadata from an MzMine 2 spectra file (.mgf)

Args: path (str): Path to spectra file (.mgf)
min_relative_intensity (float): Minimal relative intensity to keep a peak
max_relative_intensity (float): Maximal relative intensity to keep a peak
min_peaks_required (int): Minimum number of peaks to keep a spectrum
losses_from (int): minimal m/z value for losses
losses_to (int): maximal m/z value for losses
n_decimals (int): number of decimal when translating peaks/losses into words

Returns: self.document (DataFrame): A table containing spectra documents and metadata

__init__(path: str, min_relative_intensity: float = 0.01, max_relative_intensity: float = 1.0, min_peaks_required: int = 10, losses_from: int = 10, losses_to: int = 200, n_decimals: int = 2) → None

1.1 Submodules

1.1.1 memo_ms.classes module

class memo_ms.classes.**FeatureTable**(*path: str, software: str*)

Bases: `object`

Create a FeatureTable dataclass object from a feature table

Args: *path* (str): Path to a feature table file (.csv) *software* (str): One of [mzmine, xcms, msdial, memo]: the software used for feature detection.

Returns: *self.feature_table* (DataFrame): A cleaned feature quantification table

__init__(*path: str, software: str*) → `None`

export_matrix(*path, sep=','*)

Export a given matrix

Args: *path* (str): path to export *sep* (str): separator

Returns: `None`

filter(*samples_pattern, max_occurence=None*)

Filter a feature table: remove samples matching *samples_pattern* AND remove features occurring in more than *n = max_occurence* samples matched by *samples_pattern*

Args: *samples_pattern* (str): the str pattern to match in samples to filter *max_occurence* (int): maximal number of occurrence allowed in matched samples before removing a feature/word

Returns: *self.filtered_feature_table* (DataFrame): A filtered feature table

class memo_ms.classes.**MemoMatrix**

Bases: `object`

Create an empty MemoMatrix dataclass object

__init__() → `None`

export_matrix(*path, sep=','*)

Export a given matrix

Args: *path* (str): path to export *sep* (str): separator

Returns: `None`

filter(*samples_pattern, max_occurence=None*)

Filter a MEMO matrix: remove samples matching *samples_pattern* AND remove features occurring in more than *n = max_occurence* samples matched by *samples_pattern*

Args: *samples_pattern* (str): the str pattern to match in samples to filter *max_occurence* (int): maximal number of occurrence allowed in matched samples before removing a feature/word

Returns: *self.memo_matrix* (DataFrame): A filtered feature table matrix

memo_from_aligned_samples(*featuretable, spectradocuments*) → `pandas.core.frame.DataFrame`

Use a *featuretable* and a *spectradocuments* to generate a MEMO matrix. Returns a `pd.DataFrame` MEMO matrix.

Args: *featuretable* (FeatureTable): a FeatureTable dataclass object *spectradocuments* (SpectraDocuments): a SpectraDocuments dataclass object

Returns: *self.memo_matrix* (DataFrame): A MEMO matrix

```
memo_from_unaligned_samples(path_to_samples_dir, min_relative_intensity=0.01,  
                             max_relative_intensity=1.0, min_peaks_required=10, losses_from=10,  
                             losses_to=200, n_decimals=2)
```

Generate a Memo matrix from a list of individual .mgf files

Args: path_to_samples_dir (str): Path to the directory where individual .mgf files are gathered
min_relative_intensity (float): Minimal relative intensity to keep a peak
max_relative_intensity (float): Maximal relative intensity to keep a peak
min_peaks_required (int): Minimum number of peaks to keep a spectrum
losses_from (int): minimal m/z value for losses
losses_to (int): maximal m/z value for losses
n_decimals (int): number of decimal when translating peaks/losses into words

Returns: self.memo_matrix (DataFrame): A MEMO matrix

```
merge_memo(memomatrix_2, drop_not_in_common=False)
```

Merge 2 MEMO matrix

Args: memomatrix_2 (MemoContainer): MemoMatrix dataclass object containing the 2nd MEMO matrix to merge
drop_not_in_common (bool): Drop peaks/losses not in common

Returns: MemoContainer (MemoContainer): A MemoMatrix dataclass object containing the merged MEMO matrix

```
class memo_ms.classes.SpectraDocuments(path: str, min_relative_intensity: float = 0.01,  
                                       max_relative_intensity: float = 1.0, min_peaks_required: int = 10,  
                                       losses_from: int = 10, losses_to: int = 200, n_decimals: int = 2)
```

Bases: `object`

Create a SpectraDocuments dataclass object containing spectra documents and metadata from an MzMine 2 spectra file (.mgf)

Args: path (str): Path to spectra file (.mgf)
min_relative_intensity (float): Minimal relative intensity to keep a peak
max_relative_intensity (float): Maximal relative intensity to keep a peak
min_peaks_required (int): Minimum number of peaks to keep a spectrum
losses_from (int): minimal m/z value for losses
losses_to (int): maximal m/z value for losses
n_decimals (int): number of decimal when translating peaks/losses into words

Returns: self.document (DataFrame): A table containing spectra documents and metadata

```
__init__(path: str, min_relative_intensity: float = 0.01, max_relative_intensity: float = 1.0,  
         min_peaks_required: int = 10, losses_from: int = 10, losses_to: int = 200, n_decimals: int = 2)  
→ None
```

1.1.2 memo_ms.import_data module

```
memo_ms.import_data.import_memo_quant_table(path) → pandas.core.frame.DataFrame
```

Import feature quantification table memo ready

Args: path (str): Path to feature quantification table

Returns: quant_table (DataFrame): A cleaned feature quantification table

```
memo_ms.import_data.import_msdiag_quant_table(path) → pandas.core.frame.DataFrame
```

Import feature quantification table generated from MS-DIAL and clean it

Args: path (str): Path to feature quantification table

Returns: quant_table (DataFrame): A cleaned MS-DIAL feature quantification table

```
memo_ms.import_data.import_mzmine2_quant_table(path) → pandas.core.frame.DataFrame
```

Import feature quantification table generated from MzMine 2 and clean it

Args: path (str): Path to feature quantification table

Returns: quant_table (DataFrame): A cleaned MzMine2 feature quantification table

`memo_ms.import_data.import_xcms_quant_table(path) → pandas.core.frame.DataFrame`

Import feature quantification table generated from XCMS and clean it

Args: path (str): Path to feature quantification table

Returns: quant_table (DataFrame): A cleaned XCMS feature quantification table

`memo_ms.import_data.load_and_filter_from_mgf(path, min_relative_intensity, max_relative_intensity, loss_mz_from, loss_mz_to, n_required) → list`

Load and filter spectra from mgf file to prepare for MEMO matrix generation

Returns: spectrums (list of matchms.spectrum): a list of matchms.spectrum objects

1.1.3 memo_ms.visualization module

`memo_ms.visualization.plot_hca(matrix, df_metadata, filename_col, group_col, plotly_discrete_cm=['#636EFA', '#EF553B', '#00CC96', '#AB63FA', '#FFA15A', '#19D3F3', '#FF6692', '#B6E880', '#FF97FF', '#FECB52'], linkage_method='ward', linkage_metric='euclidean', norm=False, scaling=False)`

Simple HCA plot of a MEMO matrix / Feature table using matplotlib

Args: matrix (DataFrame): A Table in the MemoMatrix.memo_matrix or FeatureTable.feature_table format
df_metadata (DataFrame): Metadata of the MEMO matrix samples filename_col (str): Column name in df_metadata to match memo_matrix index group_col (str): Column name in df_metadata to use as groups for plotting plotly_discrete_cm ([type], optional): Plotly discrete colormap to use for groups. Defaults to px.colors.qualitative.Plotly. linkage_method (str, optional): Linkage method to use. Defaults to 'ward'. linkage_metric (str, optional): Linkage metric to use. Defaults to 'euclidean'. norm (bool, optional): Apply samples normalization. Defaults to False. scaling (bool, optional): Apply pareto scaling to MEMO matrix columns. Defaults to False.

Returns: None

`memo_ms.visualization.plot_heatmap(matrix, df_metadata, filename_col, group_col, plotly_discrete_cm=['#636EFA', '#EF553B', '#00CC96', '#AB63FA', '#FFA15A', '#19D3F3', '#FF6692', '#B6E880', '#FF97FF', '#FECB52'], linkage_method='ward', linkage_metric='euclidean', heatmap_metric='braycurtis', norm=False, scaling=False)`

HCA and heatmap plot of a MEMO matrix / Feature table using Plotly

Args: matrix (DataFrame): A Table in the MemoMatrix.memo_matrix or FeatureTable.feature_table format
df_metadata (DataFrame): Metadata of the MEMO matrix samples filename_col (str): Column name in df_metadata to match memo_matrix index group_col (str): Column name in df_metadata to use as groups for plotting plotly_discrete_cm ([type], optional): Plotly discrete colormap to use for groups. Defaults to px.colors.qualitative.Plotly. linkage_method (str, optional): Linkage method to use. Defaults to 'ward'. linkage_metric (str, optional): Linkage metric to use. Defaults to 'euclidean'. heatmap_metric (str, optional): Distance metric to use for heatmap. Defaults to 'braycurtis'. norm (bool, optional): Apply samples normalization. Defaults to False. scaling (bool, optional): Apply pareto scaling to MEMO matrix columns. Defaults to False.

Returns: None

`memo_ms.visualization.plot_pcoa_2d(matrix, df_metadata, filename_col, group_col, metric='braycurtis', norm=False, scaling=False, pc_to_plot=(1, 2))`

Simple 2D PCoA plot of a MEMO matrix / Feature table using Plotly

Args: matrix (DataFrame): A Table in the MemoMatrix.memo_matrix or FeatureTable.feature_table format
df_metadata (DataFrame): Metadata of the MEMO matrix samples filename_col (str): Column name in df_metadata to match memo_matrix index group_col (str): Column name in df_metadata to use as groups for plotting metric (str, optional): Distance metric to use, see <https://docs.scipy.org/doc/scipy/reference/generated/scipy.spatial.distance.pdist.html>. Defaults to 'braycurtis'. norm (bool, optional): Apply samples normalization. Defaults to False. scaling (bool, optional): Apply pareto scaling to MEMO matrix columns. Defaults to False. pc_to_plot (list of int, optional): PCs to plot. Defaults to [1,2].

Returns: None

memo_ms.visualization.plot_pcoa_3d(matrix, df_metadata, filename_col, group_col, metric='braycurtis', norm=False, scaling=False, pc_to_plot=(1, 2, 3))

Simple 2D PCoA plot of a MEMO matrix / Feature table using Plotly

Args: matrix (DataFrame): A Table in the MemoMatrix.memo_matrix or FeatureTable.feature_table format
df_metadata (DataFrame): Metadata of the MEMO matrix samples filename_col (str): Column name in df_metadata to match memo_matrix index group_col (str): Column name in df_metadata to use as groups for plotting metric (str, optional): Distance metric to use, see <https://docs.scipy.org/doc/scipy/reference/generated/scipy.spatial.distance.pdist.html>. Defaults to 'braycurtis'. norm (bool, optional): Apply samples normalization. Defaults to False. scaling (bool, optional): Apply pareto scaling to MEMO matrix columns. Defaults to False. pc_to_plot (list of int, optional): PCs to plot. Defaults to [1,2,3].

Returns: None

DESCRIPTION

MEMO is a method allowing a Retention Time (RT) agnostic alignment of metabolomics samples using the fragmentation spectra (MS2) of their constituents. The occurrence of MS2 peaks and neutral losses (to the precursor) in each sample is counted and used to generate an *MS2 fingerprint* of the sample. These fingerprints can in a second stage be aligned to compare different samples. Once obtained, different filtering (remove peaks/losses from blanks for example) and visualization techniques (MDS/PCoA, TMAP, Heatmap, ...) can be used. MEMO suits particularly well to compare chemodiverse samples, *i.e.* with a poor features overlap, or to compare samples with a strong RT shift, acquired using different LC methods or even different mass spectrometers technology (MaXis Q-ToF vs Q-Exactive).

MEMO is mainly built on [matchms](#) and [spec2vec](#) packages for handling the MS2 spectra and converting them into documents. Huge thanks to them for the amazing work done with these packages!

TO INSTALL IT:

First, make sure to have [anaconda](#) installed.

1. Create a new conda environment to avoid clashes:

```
conda create --name memo python=3.8
conda activate memo
```

2. Install with pip:

```
pip install numpy
pip install memo-ms
```

If you have an error, try installing scikit-bio from conda-forge before installing the package with pip:

```
conda install -c conda-forge scikit-bio
pip install memo-ms
```

You can clone the Github package repository to get the demo files and the tutorial!

EXAMPLES

Different examples of application and comparison to other MS/MS-based metrics are available [here](#) and the corresponding notebooks are available on [GitHub](#).

INDICES AND TABLES

- genindex
- modindex
- search

PYTHON MODULE INDEX

m

- `memo_ms`, [1](#)
- `memo_ms.classes`, [3](#)
- `memo_ms.import_data`, [4](#)
- `memo_ms.visualization`, [5](#)

Symbols

__init__() (*memo_ms.FeatureTable* method), 1
 __init__() (*memo_ms.MemoMatrix* method), 1
 __init__() (*memo_ms.SpectraDocuments* method), 2
 __init__() (*memo_ms.classes.FeatureTable* method), 3
 __init__() (*memo_ms.classes.MemoMatrix* method), 3
 __init__() (*memo_ms.classes.SpectraDocuments* method), 4

E

export_matrix() (*memo_ms.classes.FeatureTable* method), 3
 export_matrix() (*memo_ms.classes.MemoMatrix* method), 3
 export_matrix() (*memo_ms.FeatureTable* method), 1
 export_matrix() (*memo_ms.MemoMatrix* method), 1

F

FeatureTable (*class in memo_ms*), 1
 FeatureTable (*class in memo_ms.classes*), 3
 filter() (*memo_ms.classes.FeatureTable* method), 3
 filter() (*memo_ms.classes.MemoMatrix* method), 3
 filter() (*memo_ms.FeatureTable* method), 1
 filter() (*memo_ms.MemoMatrix* method), 1

I

import_memo_quant_table() (*in module memo_ms.import_data*), 4
 import_msdiag_quant_table() (*in module memo_ms.import_data*), 4
 import_mzmine2_quant_table() (*in module memo_ms.import_data*), 4
 import_xcms_quant_table() (*in module memo_ms.import_data*), 5

L

load_and_filter_from_mgf() (*in module memo_ms.import_data*), 5

M

memo_from_aligned_samples()

(*memo_ms.classes.MemoMatrix* method), 3
 memo_from_aligned_samples() (*memo_ms.MemoMatrix* method), 1
 memo_from_unaligned_samples() (*memo_ms.classes.MemoMatrix* method), 3
 memo_from_unaligned_samples() (*memo_ms.MemoMatrix* method), 2
 memo_ms
 module, 1
 memo_ms.classes
 module, 3
 memo_ms.import_data
 module, 4
 memo_ms.visualization
 module, 5
 MemoMatrix (*class in memo_ms*), 1
 MemoMatrix (*class in memo_ms.classes*), 3
 merge_memo() (*memo_ms.classes.MemoMatrix* method), 4
 merge_memo() (*memo_ms.MemoMatrix* method), 2
 module
 memo_ms, 1
 memo_ms.classes, 3
 memo_ms.import_data, 4
 memo_ms.visualization, 5

P

plot_hca() (*in module memo_ms.visualization*), 5
 plot_heatmap() (*in module memo_ms.visualization*), 5
 plot_pcoa_2d() (*in module memo_ms.visualization*), 5
 plot_pcoa_3d() (*in module memo_ms.visualization*), 6

S

SpectraDocuments (*class in memo_ms*), 2
 SpectraDocuments (*class in memo_ms.classes*), 4